



# BrainXell Motor Neuron Monoculture Protocol (v10.0)

## CONTENTS

- One vial cryopreserved human motor neurons (BrainXell BX-0100-XX-XX)
- BrainFast Motor (BrainXell # BX-2100) formerly known as Neuron Seeding Supplement
- BrainFast D4 (BrainXell # BX-2040) formerly known as Day 4 Supplement
- BrainFast SK (BrainXell # BX-2020) formerly known as Supplement K

## STORAGE

- Immediately transfer the cryovial of neurons to a liquid or vapor nitrogen storage system.
- BrainFast supplements should be stored at  $-20^{\circ}\text{C}$  (6 months) or  $-80^{\circ}\text{C}$  (18 months). Return vial to  $-20^{\circ}\text{C}$  between each time of use to maintain stability.

## ADDITIONAL MATERIALS NEEDED

- DMEM/F-12 Medium, +L-glutamine +HEPES (Thermo Fisher Scientific #11330032)
- Neurobasal Medium (Thermo Fisher Scientific #21103049)
- B-27 Supplement (Thermo Fisher Scientific #17504044)
- N-2 Supplement (Thermo Fisher Scientific #17502048)
- GlutaMAX Supplement (Thermo Fisher Scientific #35050061)
- Geltrex (Thermo Fisher Scientific #A1413201) or Cultrex (R&D Systems #3432-001-01)
- Ascorbic Acid (MilliporeSigma #A8960)
- PDL-Coated 96-Well Plates (Refer to [BrainXell Culture Plates PDL Coating Protocol v10.0](#))

## OPTIONAL MATERIALS

Not included in medium during BrainXell QC assessment

- BDNF (Thermo Fisher Scientific #450-02)
- GDNF (Thermo Fisher Scientific #450-10)
- TGF- $\beta$ 1 (Thermo Fisher Scientific #100-21C)

## PROCEDURE

### Day 0 Seeding Preparation

To seed (1) full 96-well plate you will need 2.2 – 3.3 million live cells. Additional cells may be used depending on individual research needs. Review the Certificate of Analysis (CoA) for viability and seeding information.

1. Gather the Basal Medium components according to the recipe (see Media Compositions section below).
  - a. Note that BDNF, GDNF, and TGF- $\beta$ 1 are lyophilized powders. Follow the manufacturer's instructions for reconstitution and long-term storage. We recommend creating stock solutions of 10  $\mu\text{g}/\text{mL}$  BDNF, 10  $\mu\text{g}/\text{mL}$  GDNF, and 2  $\mu\text{g}/\text{mL}$  TGF- $\beta$ 1.

- b. We recommend reconstituting ascorbic acid in DMEM/F-12 medium to create a 200 mM stock solution; See “Preparation of 200 mM Ascorbic Acid” at end of this document for important information.
2. Working in a cell culture hood (biological safety cabinet), combine all Basal Medium components in an appropriately sized sterile container. Allow the Basal Medium to equilibrate to room temperature for at least 15 minutes. Do not warm the medium in a 37°C water bath. Culture plates should also be at room temperature prior to use.
3. Prepare a pre-diluted Cultrex solution: add 495 µL of **cold** DMEM/F12 into a 55 µL aliquot of frozen Cultrex taken directly from the -80 freezer. Mix to dissolve and store at 4°C until it is time to seed the plate (Step 16).
4. Prepare one 50-mL conical tube: add 3 mL of Basal Medium.
5. Prepare one microcentrifuge tube: add 25 µL of Trypan Blue solution for cell counting.

### Day 0 Seeding the Neurons

6. Remove the cryovial of cells from nitrogen storage and place in a 37°C water bath. To minimize contamination, avoid submerging the cap. Gently move the vial within the bath to increase the rate of thawing.
7. As soon as the last of the ice melts, which will take ~75-90 seconds, remove the cryovial from the water bath. Disinfect the vial by spraying it with 70% ethanol before transferring it into the cell culture hood.
8. Using a P1000 pipette, transfer 500 µL Basal Medium from the prepared 50-mL conical tube to the cell cryovial at a rate of ~2-3 drops/sec. This process should take about 10 seconds.
9. Gently transfer all contents from the cryovial (~1 mL total) back to the same 50-mL conical tube.
10. Centrifuge cells at 1700 rpm (465xg) for 5 mins.
11. Carefully aspirate supernatant, leaving the cell pellet, and resuspend cells in 950 µL fresh Basal Medium. Mix thoroughly by gently pipetting up and down using a P1000 pipette.
12. Based on the CoA value (*Viable Cells/Vial*), resuspend the cells to a concentration of  $1.0 \times 10^6$  live cells/mL by slowly adding additional Basal Medium to the existing ~1 mL in the tube.
  - a. Example: 5.3 million Viable Cells/Vial is diluted to 5.3 mL total volume.
13. Count the cells: gently swirl the conical tube and pipette up and down 3-5 times to ensure cells are evenly suspended in medium. Transfer 25 µL of cell suspension to the microcentrifuge tube prefilled with 25 µL Trypan Blue solution from Step 5 and pipette up and down a few times to mix. Count the number of viable and dead cells using a hemocytometer. Determine the live cell concentration (live cells/mL) and viability.
14. Calculate volumes needed to make the motor neuron (MN) Seeding Suspension. A typical seeding density is 20,000 – 30,000 viable cells/100 µL/well for a 96-well plate (~62,500 – 93,750 viable cells/cm<sup>2</sup>). Recommended seeding density may vary based on lot number; refer to the CoA for lot-specific seeding information. Dilute the cells to the desired seeding concentration based on the Trypan Blue cell count.

#### Example of dilution calculations.

Actual Live Cell Concentration	Target Seeding Density	Total Seeding Volume Needed	Volume of Cell Suspension Volume	Volume of Basal Medium
Ex: $1.1 \times 10^6$ viable cells/mL (Step 13)	20,000 viable cells per 100 µL/well = 200,000 cells/mL	11 mL (Seeding 1 plate)	<b>Mix: 2 mL</b> $(\frac{200,000}{\text{mL}} \times 11\text{mL})$ $1.1 \times 10^6$ cells/mL	<b>Add: 9 mL</b> (11 mL – 2 mL)

15. In a new sterile 50-mL conical tube, mix the calculated volumes of cell suspension and Basal Medium needed to obtain a final volume of 11 mL MN Seeding Suspension.
16. Add appropriate volumes of BrainFast Motor (1:1000) and pre-diluted Cultrex from Step 3 (1:100) to the MN Seeding Suspension per the Seeding Medium recipe.
  - a. Example: add 11  $\mu$ L of BrainFast Motor and 110  $\mu$ L pre-diluted Cultrex solution to 11 mL of MN Seeding Suspension.
17. Mix completely and transfer 100  $\mu$ L/well (Ex: 20,000 cells/well) of the final Seeding Medium into a PDL-coated 96-well plate using a multi-channel pipette or liquid handler. Do not move or agitate the plate throughout the duration of the seeding process as this may lead to uneven attachment.
18. After seeding, let the plate rest for 10 minutes before moving it to allow the cells to settle to the bottom of the well. After 10 minutes, very gently transfer the plate to a humidified incubator at 37°C with 5% CO<sub>2</sub>. The day of cell plating (today) is designated as Day 0.

*\*Note: The entire thawing and plating process should not exceed 1 hour. Post-thaw viability and overall cell health could be severely impacted and lead to an unsuccessful culture if the whole process is too long.*

### Day 1 Medium Replacement

19. On Day 1 (24 hours after seeding), prepare fresh Day 1 Medium.
20. Carefully remove all 100  $\mu$ L medium/well and gently replace with 100  $\mu$ L/well Day 1 Medium. Complete one row or column at a time to ensure the wells do not dry out during the medium replacement process.

### Day 4 Medium Addition

21. On Day 4 (96 hours after seeding), prepare fresh Day 4 Medium.
  - a. Note: If only culturing for 7 days the addition of BrainFast SK on Day 4 is optional for motor neurons. If culturing longer than 7 days, the addition of BrainFast SK on Day 4 is recommended.
22. Gently add 100  $\mu$ L/well of Day 4 medium to the entire plate for a total of 200  $\mu$ L/well.

### Day 7 and Onward Medium Changes

23. Change half the medium weekly (on Day 7, 14, 21, etc.) using Basal Medium (made in Step 2).
  - a. Note: Addition of low concentration SK (0.1X-0.5X) in medium may be helpful for long-duration cultures; please contact [support@brainxell.com](mailto:support@brainxell.com) for further assistance.
  - b. Gently remove 100  $\mu$ L/well and slowly add 100  $\mu$ L/well Basal Medium to the entire plate.

The neurons mature rapidly and can be maintained viable and adherent in culture under the above conditions for at least 3 weeks post-seeding.

## Motor Neuron Media Compositions

Step 2	Component	Stock Conc.	Final Conc.	1 Plate Volume	2 Plate Volume	5 Plate Volume	
Basal Medium	1	DMEM/F12 Medium	1X	0.5X	24 mL	48 mL	120 mL
	2	Neurobasal Medium	1X	0.5X	24 mL	48 mL	120 mL
	3	B-27 Supplement	50X	1X	1 mL	2 mL	5 mL
	4	N-2 Supplement	100X	1X	500 µL	1 mL	2.5 mL
	5	GlutaMAX	200 mM	0.5 mM	125 µL	250 µL	625 µL
	6	Ascorbic Acid Solution (Step1b)	200 mM	0.2 mM	50 µL	100 µL	250 µL
Basal Medium stock can be stored for up to 3 weeks at 4°C.							

Note: Addition of Growth Factors is optional. If added, recommended final medium concentrations are 10 ng/mL BDNF, 10 ng/mL GDNF, and 1 ng/mL TGF-β1.

Step 16	Component	Stock Conc.	Final Conc.	1 Plate Volume	2 Plate Volume	5 Plate Volume	
Seeding Medium	1	<b>MN Seeding Suspension</b> (Step 15)	<b>1X</b>	<b>1X</b>	<b>11 mL</b>	<b>22 mL</b>	<b>55 mL</b>
	2	BrainFast Motor	1000X	1X	11 µL	22 µL	55 µL
	3	Pre-diluted Cultrex*	1.5mg/mL	15 µg/mL	110 µL	220 µL	550 µL

\* Cultrex must be pre-diluted according to step 3 in the protocol and added to the final seeding suspension right before seeding the plate.

Step 19	Component	Stock Conc.	Final Conc.	1 Plate Volume	2 Plate Volume	5 Plate Volume	
Day 1 Medium	1	<b>Basal Medium</b>	<b>1X</b>	<b>1X</b>	<b>11 mL</b>	<b>22 mL</b>	<b>55 mL</b>
	2	BrainFast Motor	1000X	1X	11 µL	22 µL	55 µL
	3	BrainFast SK	1000X	1X	11 µL	22 µL	55 µL

Step 21	Component	Stock Conc.	Final Conc.	1 Plate Volume	2 Plate Volume	5 Plate Volume	
Day 4 Medium	1	<b>Basal Medium</b>	<b>1X</b>	<b>1X</b>	<b>11 mL</b>	<b>22 mL</b>	<b>55 mL</b>
	2	BrainFast D4	1000X	1X	11 µL	22 µL	55 µL
	3	BrainFast SK*	1000X	1X	11 µL	22 µL	55 µL

\* BrainFast SK addition is optional in Day 4 Medium for cultures that will be grown 7 days or less.

## Preparation of 200 mM Ascorbic Acid Solution

### Additional Sterile Materials Required

50-mL Conical Tubes, 0.22 micron Filter, Amber Microcentrifuge Tubes

### Critical Information

- Ascorbic Acid must be added *into* the DMEM/F12 (not vice versa) to ensure it dissolves.
- The rate of Ascorbic Acid addition to DMEM/F12 is critical. If added too rapidly, the chemical will fail to go into solution.
- Only mix by inversion. Undissolved Ascorbic Acid will adhere to pipette tips and other implements.
- Ascorbic Acid is pH, temperature, and light sensitive.

### Steps

1. Aliquot 17.2 mL of DMEM/F12 into a 50-ml tube.
2. Weigh 1.0 gram of Ascorbic Acid using a separate vessel.
3. **Critical Step:** Gradually transfer Ascorbic Acid **into** the tube of DMEM/F12 (no more than 0.2 gram at a time, waiting 10 seconds between each addition).
4. Secure the cap on the tube and invert to mix.
5. Allow the Ascorbic Acid to dissolve into the solution:
  - a. Place in a 37°C water bath for 30 minutes, inverting the tube every 10 minutes to mix.
  - b. After 30 minutes check to see if Ascorbic Acid is fully dissolved. If not, invert tube to mix and place capped tube at 4°C overnight.
  - c. After overnight storage, invert to mix and confirm Ascorbic Acid is completely dissolved.
6. Filter the fully dissolved Ascorbic Acid Solution using a Steriflip filter or equivalent 0.22 µm sterile filter.
7. Aliquot into amber 1.5mL Eppendorf tubes; recommended volumes of 110µl or 550µl.
8. Maintain frozen aliquots at -20°C for up to 6 months. Aliquots are single use. Do not refreeze.